

# Overdenture VulkanLoc®



## VulkanLoc® Abutment (Grade 5 Titanium + TiN Coating)


30 Ncm (max.)

|    |             |             |             |             |             |             |             |
|----|-------------|-------------|-------------|-------------|-------------|-------------|-------------|
|    |             |             |             |             |             |             |             |
|    | ±0.7mm      | 2mm         | 3mm         | 4mm         | 5mm         | 6mm         | 1.5mm       |
|    | h1.0 mm     | h2.0 mm     | h3.0 mm     | h4.0 mm     | h5.0 mm     | h6.0 mm     | Soft Tissue |
| NP | VEXT11N1000 | VEXT11N2000 | VEXT11N3000 | VEXT11N4000 | VEXT11N5000 | -           |             |
| RP | VEXT11R1000 | VEXT11R2000 | VEXT11R3000 | VEXT11R4000 | VEXT11R5000 | VEXT11R6000 |             |

## Retentions VulkanLoc®

|   |   |
|---|---|
| <ul style="list-style-type: none"> <li>Blue Ref: L0402 X4 (1.5 lbs)</li> <li>Pink Ref: L0403 X4 (3.0 lbs)</li> <li>Transparent Ref: L0404 X4 (5.0 lbs)</li> </ul> | <ul style="list-style-type: none"> <li>Black Ref: L0400 X4</li> <li>Red Ref: L0401 X4 (1.0 lbs)</li> <li>Orange Ref: L0405 X4 (2.0 lbs)</li> <li>Green Ref: L0407 X4 (4.0 lbs)</li> </ul> |
|---|---|

## VulkanLoc® Processing Package



## Housing with Black Retainer



## VulkanLoc® Analog



## Mounter for VulkanLoc®



## Spacer Ring



## Smart Tool VulkanLoc®



Multi Functional Tool for VulkanLoc® System

## VulkanLoc® Driver



## Imp. Coping for VulkanLoc®

