


Vulkan Tissue Care

*Only for single prostheses



Straight Tissue Care Abutment + Screw



Tissue Care Screw



Healing

Healing Cap for Tissue Care



Impression

Impression Coping for Tissue Care



Analog for Tissue Care



Screwed prosthesis

Castable Tissue Care



Cobalt-Chrome Castable Tissue Care



Temporary Abutment Tissue Care



Clinical Screw Tissue Care



Angled Prosthesis

Angled Cobalt-Chrome



- 1 VCHA170001
- 2 VCHA300001
- 3 TCS0500-A

Angled Screw

