



CONICAL CONNECTION



PRÓTESIS ATORNILLADA

Unitarios
(Referencias Antirrotatorias)

¿Usa Transepiteliales?

Tissue Care
H:1,5mm. Ref: VCON14N0015-X
H:2,5mm. Ref: VCON14N0025-X
H:3,5mm. Ref: VCON14N0035-X
H:4,5mm. Ref: VCON14N0045-X

Prot. Recta
Prot. Angulada

Coping Impresión.
Ref: TCS0201

Calcinable. Ref: TCS0400
Cr-Co. Ref: TCS0500
Provisional. Ref: TCS0600

TORNILLO CLÍNICO Y LLAVE PROTÉSICA

TITANIO Ref: TCS0957
DLC Ref: TCS0957T

LLAVE UNIGRIP
Ref: VSDUG-1-R Ref: VSDUG-3
Ref: VSDUG-2-R Ref: VSDUG-4

TITANIO TLB Ref: TCS0900-TLB

LLAVE TETRALOBULAR
Ref: VSDTLB-1 Ref: VSDTLB-3
Ref: VSDTLB-2 Ref: VSDTLB-4

TITANIO Ref: VCON09N07
DLC Ref: VCON09N07T

LLAVE UNIGRIP
Ref: VSDUG-1-R Ref: VSDUG-3
Ref: VSDUG-2-R Ref: VSDUG-4

TITANIO TLB Ref: VCON09N00-TLB

LLAVE TETRALOBULAR
Ref: VSDTLB-1 Ref: VSDTLB-3
Ref: VSDTLB-2 Ref: VSDTLB-4

Calcinable. Ref: MU0402
Cr-Co. Ref: MU0502
Provisional. Refs: MU0602 (Titanio) | MU0602P (Peek)

TITANIO Ref: MU0905
DLC Ref: MU0905T

LLAVE UNIGRIP
Ref: VSDUG-1-R Ref: VSDUG-3
Ref: VSDUG-2-R Ref: VSDUG-4

TITANIO TLB Ref: MU0900-TLB

LLAVE TETRALOBULAR
Ref: VSDTLB-1 Ref: VSDTLB-3
Ref: VSDTLB-2 Ref: VSDTLB-4

SOLUCIÓN ANGULABLE
Cr-Co. Angulado Ref: MU0504

Coping Impresión.
Cubeta Abierta: Ref: MU0201 | Cubeta Cerrada: Ref: MU0202

Múltiples
(Referencias Rotatorias)

¿Usa Transepiteliales?

Prot. Recta

Prot. Angulada

Multi-Use Recto
H:1,5mm. Ref: VCON10N1500
H:2,5mm. Ref: VCON10N2500
H:3,5mm. Ref: VCON10N3500
H:4,5mm. Ref: VCON10N4500
H:5,5mm. Ref: VCON10N5500
H:6,5mm. Ref: VCON10N6500

Calcinable. Ref: MU0402
Cr-Co. Ref: MU0502
Provisional. Refs: MU0602 (Titanio) | MU0602P (Peek)

TITANIO Ref: MU0905
DLC Ref: MU0905T

LLAVE UNIGRIP
Ref: VSDUG-1-R Ref: VSDUG-3
Ref: VSDUG-2-R Ref: VSDUG-4

TITANIO TLB Ref: MU0900-TLB

LLAVE TETRALOBULAR
Ref: VSDTLB-1 Ref: VSDTLB-3
Ref: VSDTLB-2 Ref: VSDTLB-4

SOLUCIÓN ANGULABLE
Cr-Co. Angulado Ref: MU0504

Prot. Recta

Prot. Angulada

Multi-Use Angulado
17° H:2,5mm. Ref: VCON10N2517
17° H:3,5mm. Ref: VCON10N3517
30° H:3,5mm. Ref: VCON10N3530
30° H:4,5mm. Ref: VCON10N4530

Calcinable. Ref: MU0402
Cr-Co. Ref: MU0502
Provisional. Refs: MU0602 (Titanio) | MU0602P (Peek)

TITANIO Ref: MU0905
DLC Ref: MU0905T

LLAVE UNIGRIP
Ref: VSDUG-1-R Ref: VSDUG-3
Ref: VSDUG-2-R Ref: VSDUG-4

TITANIO TLB Ref: MU0900-TLB

LLAVE TETRALOBULAR
Ref: VSDTLB-1 Ref: VSDTLB-3
Ref: VSDTLB-2 Ref: VSDTLB-4

SOLUCIÓN ANGULABLE
Cr-Co. Angulado Ref: MU0504

Prot. Recta

Prot. Angulada

Calcinable. Ref: VCON04N0002
Cr-Co. Ref: VCON05N0002
Provisional. Ref: VCON06N0002

TITANIO Ref: VCON09N07
DLC Ref: VCON09N07T

LLAVE UNIGRIP
Ref: VSDUG-1-R Ref: VSDUG-3
Ref: VSDUG-2-R Ref: VSDUG-4

TITANIO TLB Ref: VCON09N00-TLB

LLAVE TETRALOBULAR
Ref: VSDTLB-1 Ref: VSDTLB-3
Ref: VSDTLB-2 Ref: VSDTLB-4

17° Ref: VCHA170001
30° Ref: VCHA170001

PRÓTESIS CEMENTADA

Prótesis Unitaria
Prótesis Múltiple

Prot. Recta

Prot. Angulada

H:1,5mm. Ref: VCON07N1501
H:2,5mm. Ref: VCON07N2501
H:3,5mm. Ref: VCON07N3501
H:4,5mm. Ref: VCON07N4501

H:1,5mm. Ref: VCON07N1515
H:2,5mm. Ref: VCON07N2515
H:3,5mm. Ref: VCON07N3515

H:1,5mm. Ref: VCON07N1525
H:2,5mm. Ref: VCON07N2525
H:3,5mm. Ref: VCON07N3525

TITANIO Ref: VCON09N07
DLC Ref: VCON09N07T

LLAVE UNIGRIP
Ref: VSDUG-1-R Ref: VSDUG-3
Ref: VSDUG-2-R Ref: VSDUG-4

PRÓTESIS REMOVIBLE

VulkanLoc
Prótesis para barra

Locator Recto

H:1,0mm. Ref: VCON11N1000
H:2,0mm. Ref: VCON11N2000
H:3,0mm. Ref: VCON11N3000

H:4,0mm. Ref: VCON11N4000
H:5,0mm. Ref: VCON11N5000
H:6,0mm. Ref: VCON11N6000

Mirar Prótesis Atornillada Múltiple

Unitario
(Referencias Antirrotatorias)

¿Usa Transepiteliales?

Tissue Care
Todos lo Pilares
Tissue Care

Scan Body Tissue Care Ref: TCS1301
Interfase Ang. Tissue Care H:0,5mm Ref: TCS080501
Análogo Digital Tissue Care Ref: TCS0310

Prot. Recta
Prot. Angulada

TITANIO Ref: TCS0957
DLC Ref: TCS0957T

LLAVE UNIGRIP
Ref: VSDUG-1-R Ref: VSDUG-3
Ref: VSDUG-2-R Ref: VSDUG-4

TITANIO TLB Ref: TCS0900-TLB

LLAVE TETRALOBULAR
Ref: VSDTLB-1 Ref: VSDTLB-3
Ref: VSDTLB-2 Ref: VSDTLB-4

Scan Body Directo Implante Ref: VCON13N0002
Ref: VCON09N07

Interfase Angulable
H:1,5mm Ref: VCON08N1501
H:2,5mm Ref: VCON08N2501
H:3,5mm Ref: VCON08N3501

Prot. Recta
Prot. Angulada

TITANIO Ref: VCON09N07
DLC Ref: VCON09N07T

LLAVE UNIGRIP
Ref: VSDUG-1-R Ref: VSDUG-3
Ref: VSDUG-2-R Ref: VSDUG-4

TITANIO TLB Ref: VCON09N00-TLB

LLAVE TETRALOBULAR
Ref: VSDTLB-1 Ref: VSDTLB-3
Ref: VSDTLB-2 Ref: VSDTLB-4

Múltiple
(Referencias Rotatorias)

¿Usa Transepiteliales?

Multi Use
Todos lo Pilares
Multi Use

Scan Body Multi Use Ref: VMU13R02
Interfase Ang. Multi Use H:0,5mm Ref: MU080502
Análogo Digital Multi Use Ref: MU0310

Prot. Recta
Prot. Angulada

TITANIO Ref: MU0905
DLC Ref: MU0905T

LLAVE UNIGRIP
Ref: VSDUG-1-R Ref: VSDUG-3
Ref: VSDUG-2-R Ref: VSDUG-4

TITANIO TLB Ref: MU0900-TLB

LLAVE TETRALOBULAR
Ref: VSDTLB-1 Ref: VSDTLB-3
Ref: VSDTLB-2 Ref: VSDTLB-4

Scan Body Directo Implante Ref: VCON13N0002
Ref: VCON09N07

Interfase Angulable
H:1,5mm Ref: VCON08N1502
H:2,5mm Ref: VCON08N2502
H:3,5mm Ref: VCON08N3502

Prot. Recta
Prot. Angulada

TITANIO Ref: VCON09N07
DLC Ref: VCON09N07T

LLAVE UNIGRIP
Ref: VSDUG-1-R Ref: VSDUG-3
Ref: VSDUG-2-R Ref: VSDUG-4

TITANIO TLB Ref: VCON09N00-TLB

LLAVE TETRALOBULAR
Ref: VSDTLB-1 Ref: VSDTLB-3
Ref: VSDTLB-2 Ref: VSDTLB-4